



CARDIAC TROPONIN PANEL # MDXPCT- 004

0.5 ml Matched Serum, K2 EDTA Plasma, and Li Heparin Plasma

ID #	Age	Gender	Ethnicity	Race	Smoker/Non-Smoker	Cause of elevated Troponin	Symptoms	Date of Event	Admitted AM/PM
7006	61	Male	Hispanic	Caucasian	Non-Smoker	Heart Attack	Chest pain, left arm numbness	4/1/2019	PM

Sample ID	Draw	Time Point	BP	Draw Date	Troponin	Myoglobin	Matrix
A11703195	1	12-18 HRS	160/100	4/2/2019	40.0 ng/mL	123.5 ng/mL	SERUM
A11703196	1	12-18 HRS	160/100	4/2/2019	Not Tested	Not Tested	K2 EDTA
A11703197	1	12-18 HRS	160/100	4/2/2019	Not Tested	Not Tested	LIH
A11703201	2	24 HRS POST	150/90	4/3/2019	36.4 ng/ML	27.0 ng/mL	SERUM
A11703202	2	24 HRS POST	150/90	4/3/2019	Not Tested	Not Tested	K2 EDTA
A11703203	2	24 HRS POST	150/90	4/3/2019	Not Tested	Not Tested	LIH
A11703207	3	48 HRS POST	140/80	4/4/2019	30.0 ng/ML	13.7 ng/mL	SERUM
A11703208	3	48 HRS POST	140/80	4/4/2019	Not Tested	Not Tested	K2 EDTA
A11703209	3	48 HRS POST	140/80	4/4/2019	Not Tested	Not Tested	LIH
A11703213	4	72 HRS POST	138/82	4/5/2019	24.2 ng/ML	4.4 ng/mL	SERUM
A11703214	4	72 HRS POST	138/82	4/5/2019	Not Tested	Not Tested	K2 EDTA
A11703215	4	72 HRS POST	138/82	4/5/2019	Not Tested	Not Tested	LIH

Troponin Testing was performed on Tosoh AIA-360/Myoglobin Testing was performed on RayBio Human Myoglobin ELISA Kit

CAUTION: HUMAN BLOOD

All materials listed above are considered to be potentially infectious and should be handled with universal precautions.