



CARDIAC TROPONIN PANEL # MDXPCT- 010
0.5 ml Matched Serum, K2 EDTA Plasma, and Li Heparin Plasma

ID #	Age	Gender	Ethnicity	Race	Smoker/Non-Smoker	Cause of elevated Troponin	Symptoms	Date of Event	Admitted AM/PM
7014	57	Male	Hispanic	Caucasian	Smoker	N/A	Chest pain, left arm pain	4/17/2019	PM

Sample ID	Draw	Time Point	BP	Draw Date	Troponin	Myoglobin	Matrix
A11703288	1	12-18 HRS	178/106	4/18/2019	38.2 ng/mL	20.5 ng/mL	SERUM
A11703289	1	12-18 HRS	178/106	4/18/2019	Not Tested	Not Tested	K2 EDTA
A11703290	1	12-18 HRS	178/106	4/18/2019	Not Tested	Not Tested	LIH
A11703294	2	24 HRS POST	168/100	4/19/2019	30.0 ng/mL	2.9 ng/mL	SERUM
A11703295	2	24 HRS POST	168/100	4/19/2019	Not Tested	Not Tested	K2 EDTA
A11703296	2	24 HRS POST	168/100	4/19/2019	Not Tested	Not Tested	LIH
A11703309	3	48 HRS POST	150/90	4/20/2019	26.0 ng/mL	2.2 ng/mL	SERUM
A11703310	3	48 HRS POST	150/90	4/20/2019	Not Tested	Not Tested	K2 EDTA
A11703311	3	48 HRS POST	150/90	4/20/2019	Not Tested	Not Tested	LIH
A11703324	4	72 HRS POST	140/90	4/21/2019	24.2 ng/mL	1.4 ng/mL	SERUM
A11703325	4	72 HRS POST	140/90	4/21/2019	Not Tested	Not Tested	K2 EDTA
A11703326	4	72 HRS POST	140/90	4/21/2019	Not Tested	Not Tested	LIH

Troponin Testing was performed on Tosoh AIA-360/Myoglobin Testing was performed on RayBio Human Myoglobin ELISA Kit

CAUTION: HUMAN BLOOD

All materials listed above are considered to be potentially infectious and should be handled with universal precautions.